


Trixylyl phosphate (TXP)

BMF 79 - Trixylyl phosphate (TXP)

A SUBSTANCE OF VERY HIGH CONCERN (SVHC) BECAUSE OF ITS CMR* PROPERTIES



Toxicity:

Trixylyl phosphate (CAS 25155-23-1) is generally used as a flame retardant in hydraulic fluids and as a plasticizer.

According to Regulation (EC) No. 1272/2008 (Annex VI, Part 3) trixylyl phosphate meets the criteria for classification as a reproductive toxin (Rep. 1B (H360F: "May damage fertility") in accordance with Article 57(c) of REACH (Regulation (EC) 1907/2006).

Background:

Trixylyl phosphate (TXP) is a UVCB substance (substance of Unknown or Variable composition, Complex reaction products, or Biological materials), thus contains over 50 different constituents. It has no additives and is produced through the reaction of phosphorus oxytrichloride and xylenols. The xylenols (dimethylphenols) are present in a distillation fraction of naturally occurring coal tar derivatives, which also contains different ethyl phenols. The reaction of the different xylenols and ethylphenols results in alkylated triphenyl phosphates with a high amount of possible isomers that cannot be easily analysed according to ECHA (European chemical agency).

Therefore, the exact composition is unknown. In studies that form the basis for harmonized classification and labelling, the substance has been tested as such, but the individual constituents are not the basis for the classification.

Nobile (1980)¹ investigated the composition of two TXP products, and found the xylyl phosphate isomers 2,5-, 2,3-, 3,5-, 2,4- and 3,4- in decreasing order of abundance. The 2,6-isomer was not present. Other components identified included 4-ethylphenol, p-cresol, phenol, and trimethyl phenol.

*CMR: Carcinogenic, mutagenic, or toxic for reproduction.



BMF 79 - Trixylyl phosphate (TXP)

Preparation of reference materials:

Two samples of commercial TXP were tested and shown to have similar GC-MS and LC-UV chromatograms.

The commercial trixylyl phosphate mixture was found to contain not only TXP, but a series of compounds with higher and lower molecular weights. The main task was to purify a sample containing the mass m/z 410 from the mixture. Focusing on the molecular weight distinction, the TXP mixture was successfully separated into five fractions, the one with the expected mass m/z 410, one with a higher mass containing one more methyl group, and three with lower masses, lacking 1-3 methyl group compared to that of TXP.

The LC-UV chromatograms of the mixture and the separated fractions are shown in *Fig. 1a and 1b*.

The GC-MS chromatogram of the technical mixture before separation is shown in *Fig. 2*. The GC-MS analysis and the mass spectra of the main peaks of each of the fractions are shown in *Fig. 3*.

Available standards:

These individual TXP components are now available as reference materials from Chiron.**

Chiron No.	Description	m/z	Abbreviation
Separated fraction:			
11214.21-100-IO *	Trixylylphosphate, purified Ion m/z 368 (trimethyl)	m/z 368	TXP-3xMethyl
11215.22-100-IO	Trixylylphosphate, purified Ion m/z 382 (tetramethyl)	m/z 382	TXP-2xMethyl
11216.23-100-IO	Trixylylphosphate, purified Ion m/z 396 (pentamethyl)	m/z 396	TXP-1xMethyl
10279.24-100-IO	Trixylylphosphate, purified Ion m/z 410 (hexamethyl)	m/z 410	TXP purified
11217.25-100-IO	Trixylylphosphate, purified Ion m/z 424 (heptamethyl)	m/z 424	TXP+1xMethyl
Technical mixture:			
10265.24-100-IO	Trixylylphosphate, technical mixture		TXP
10265.24-K-IO	Trixylylphosphate, technical mixture		TXP
10265.24-50MG	Trixylylphosphate, technical mixture		TXP
Individual compounds:			
2135.21-100-IO	Tris(2-methylphenyl) phosphate (ortho)	m/z 368	2-TTP
2134.21-100-IO	Tris(3-methylphenyl) phosphate (meta)	m/z 368	3-TTP
2136.21-100-IO	Tris(4-methylphenyl) phosphate (para)	m/z 368	4-TTP
2137.21-100-IO	Tris(methylphenyl) phosphate, mix of isomers	m/z 368	TTP
11143.24-100-IO	Tris(2,3-dimethylphenyl) phosphate	m/z 410	2,3-TXP
11144.24-100-IO	Tris(2,4-dimethylphenyl) phosphate	m/z 410	2,4-TXP
11145.24-100-IO	Tris(2,5-dimethylphenyl) phosphate	m/z 410	2,5-TXP
11146.24-100-IO	Tris(2,6-dimethylphenyl) phosphate	m/z 410	2,6-TXP
11147.24-100-IO	Tris(3,4-dimethylphenyl) phosphate	m/z 410	3,4-TXP
11148.24-100-IO	Tris(3,5-dimethylphenyl) phosphate	m/z 410	3,5-TXP

Ref: (1) http://echa.europa.eu/documents/10162/19377756/supdoc_trixylyl_phosphate_en.pdf

*100-IO = 100µg/mL in isoctane, K-IO = 1000µg/mL in isoctane

**For products, pack sizes and presentations not listed, please contact sales@chiron.no



BMF 79 - Trixylyl phosphate (TXP)

Figure 1a: LC-UV chromatogram of technical trixylyl phosphate

Instrument: Agilent 1290 Infinity, UHPLC

Column: Restek, Raptor Biphenyl, 2,7 μm , 100 x 2,1 mm.

Mobile phase: Acidic gradient, Water (0,1%, v/v Formic acid + 2 mM Ammonium formate) + Methanol (0,1%, v/v Formic acid + 2 mM Ammonium formate)

Detector: DAD (λ 210 nm)





Figure 1b: LC-UV after separation of 5 fractions of the technical trixylyl phosphate





BMF 79 - Trixylyl phosphate (TXP)

Figure 2: GC-MS of technical trixylyl phosphate






Figure 3: GC-chromatograms and mass spectra of 5 trixylyl phosphate fractions (of peak *)



11214.21

TXP-3xMethyl (m/z 368)



11215.22



TXP-2xMethyl (m/z 382)



BMF 79 - Trixylyl phosphate (TXP)



11216.23

TXP-1xMethyl (m/z 398)





10279.24

TXP purified (m/z 410)



11217.25

TXP+1 xMethyl (m/z 428)





BMF 79 - Trixylyl phosphate (TXP)



*For ordering and information about prices and delivery in your country, please contact your **local distributor**:*